

Moore Balance Brace (MBB) Rx Form

Account # _____ P.O. # _____ Patient Name: _____
Account Name: _____ Date of Birth: ___/___/___ Sex: M F
E-Mail: _____ Height: _____ Weight: _____ Shoe Size: _____
Practitioner: _____ SHOES ENCLOSED
Phone: _____ Fax: _____ Primary Activity for Orthotic Use: _____
Address: _____
City, ST/Province: _____ Diagnosis: _____
Zip/Postal Code: _____

MOORE BALANCE BRACE (MBB)



MBB Standard



MBB Articulated

MBB Closure Options

- Velcro® Straps
 Lace-up

Note

- MBB are prescribed bilaterally
- A pair of flat PPT® insoles are included with every MBB order and should be worn with the AFO's

Required Measurements:

Width of forefoot at metatarsal heads _____ in Circumference above ankle _____ in
Width of ankle joint at widest point _____ in Circumference at 3" above ankle _____ in

All braces ship with a free STS Casting Sock. If no size is selected, a large sock will be sent.

Small Medium Large X-Large

Lab Cast Correction Fee:

Correction to the AFO after fabrication will incur charge. To avoid additional charges, casts must be taken with the ankle at 90 degrees neutral and the foot in subtalar neutral.

MBB Considerations:

- The MBB Fall Risk Assessment document should be performed to determine fall risk
- Medicare requires there to be a prescription or an order for the MBB in the patient chart as well as documentation of Medical Necessity, Certificate of Receipt and dispensing chart notes
- A physical/occupational therapy referral is strongly recommended for further assessment and treatment for imbalance
- Mark the cast carefully for any abnormal prominences present including malleoli
- Manually reduce any flexible deformity in casting to bring the cast to perpendicular to ensure the best fit